

**1. Part Number:** 8045P-1000

**2. Photograph:**



**3. Product Information:**

|                                 |  |
|---------------------------------|--|
| Description:                    | SteriWare PharmaScoop® PP 1000ml   |
| Use:                            | Single use, disposable scoop   |
| Material of Construction:       | Polypropylene – rigid, opaque <ul style="list-style-type: none"> <li>• FDA acceptable CFR 177.1520</li> <li>• EU Regulations 10/2011 compliant</li> <li>• EC regulation 1935/2004 compliant</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:         | Single piece, injection moulded  |
| Nominal Scoop Volume:           | 1000ml   |
| Overall Length:                 | 324mm  |
| Handle Length:                  | 149mm  |
| Maximum Width:                  | 150mm  |
| Nominal Weight:                 | 138g   |
| Moulding Environment:           | Class 100,000 Cleanroom  |
| Individually Bagged?            | Yes (heat sealed PE/PET bag)   |
| Method of Sterilisation:        | N/A  |
| Number of Scoops per Box:       | 100 (supplied as 4 boxes of 25 scoops)   |
| Recommended Storage Conditions: | Dry and ambient temperature  |
| Shelf Life:                     | 5 years from the date of manufacture   |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Dec-19